

Baby Hippy Trainer

Directions for Use



Laerdal[™]
helping save lives

Table of Contents

Instructions for Use	4
Care and Maintenance	5

Baby Hippy (375-35001) is a lifelike model of the lower torso and limbs of a female newborn. It is specifically designed for training professionals in the practice of the Ortolani Jerk-Sign and the Barlow Maneuver.

The manufacturing quality of this trainer should provide many sessions of training when reasonable care and maintenance is practiced.

Item Included:

- (1) Baby Hippy
- (1) Directions for Use
- (1) Laerdal Global Warranty Booklet

Instructions for Use:

1. Left hip is dislocated. Ortolani Jerk-Sign will reduce the femoral head into the hip socket.
2. The right hip is “lax.” Can be dislocated **ONLY** by performing the Barlow Maneuver.

WARNING: Baby Hippy is not designed to test full abduction of the hip in flexion. Attempts to fully abduct the hip will cause tearing of the simulator’s skin. Such tearing from misuse will void the warranty.

Care and Maintenance:

1. Clean with mild soap and water; do not submerge trainer in cleaning fluids or water.
2. Use only on clean surface. Avoid felt tipped markers, ink pens, acetone, iodine or other staining products and avoid placing the trainer on newsprint or inked lines of any kind.
3. To ensure longevity, trainer should be cleaned after each training session and a general inspection should be conducted regularly.
4. Modules and all other parts should be drained and air-dried thoroughly before storage and disinfected when needed.
5. Articulating parts will benefit from a light application of talcum powder prior to training sessions.
6. Store properly between teaching sessions.



Cautions and Warnings

This product contains **Natural Rubber** latex which may cause allergic reactions when in contact with humans.

© Laerdal Medical, Gatesville TX. All rights reserved
1001654 rev B
August 30, 2006