

# Neonatal Intubation Trainer

Directions for Use



**Table of Contents**

Introduction.....2  
 Laerdal Recommends .....2  
 Items Included.....2  
 Skills Taught.....2  
 Instructions for Use.....2  
 Care and Maintenance.....2  
 Replacement Parts .....2

**Introduction**

The Laerdal™ Neonatal Intubation Trainer (250-00101) is a lifelike manikin head which realistically simulates that of a newborn baby. It is specifically designed for training professionals in the practice of oral and nasal intubation.

The manufacturing quality and ease of replacing individual parts of this trainer should provide many years of training when reasonable care and maintenance is practiced.

**Laerdal Recommends**

- Uncuffed endotracheal tube (internal diameter 2.5mm)
- Manikin Lubricant or liquid soap
- Laryngoscope handle
- Laryngoscope blade #0 or #1



**Caution: Latex**

This product contains **Natural Rubber** latex which may cause allergic reactions when in contact with humans.

**Items Included**

- (1) Laerdal™ Neonatal Intubation Trainer
- (1) Manikin Lubricant
- (1) Laerdal Global Warranty Booklet
- (1) Directions for Use

**Skills Taught**

- Oral and nasal intubation
- Bag-Valve Mask ventilation
- Correct Tube placement

**Instructions for Use**

For intubation, the following equipment is needed:

- Uncuffed endotracheal tube (internal diameter 2.5mm)
- Lubricant ( or liquid soap)
- Laryngoscope handle
- Laryngoscope blade #0 or #1

1. Spray all tubes with Manikin Lubricant prior to airway placement.
2. Spray the inside of each nostril with Manikin Lubricant prior to nasal intubation.

**Care and Maintenance**

1. Clean with mild soap and water; do not submerge manikin or parts in cleaning fluids or water.
2. Use only on clean surface. Avoid felt tipped markers, ink pens, acetone, iodine or other staining products and avoid placing the manikin on newsprint or inked lines of any kind.
3. To ensure longevity, each manikin and module should be cleaned after each training session and a general inspection should be conducted regularly.
4. Modules and all other parts should be drained and air-dried thoroughly before storage and disinfected when needed.
5. Articulating parts will benefit from a light application of talcum powder prior to training sessions.
6. Store properly between teaching sessions.

**Replacement Parts**

Product Number:	Description:
250-00350.....	Lung Set
250-00450.....	Head Assembly
250-00550.....	Base/Cover Assembly
300-00650.....	Manikin Lubricant

Please contact Customer Service Representatives for more information on Replacement Parts in other skin tones.



This product is in compliance with the essential requirements of Council Directive 89 / 336 / EEC; EMC Directive.

**Notes:**



**Laerdal**

*helping save lives*