

Saver One T[®]

AED TRAINER new system



 **ami.italia**
medical products

*the ultimate preparation tool
to empower your students
to save lives*



Saver One T

Realistic and affordable AED Trainer which simulates the Saver One AED in appearance, operation, and voice prompts. Upgraded with new AHA/ERC instructions.

Extremely safe (no high-voltage), it provides realistic and versatile training for all skill levels with easy-to-follow lighting directions. Designed to meet the needs of all trainers, for teaching both first time lay users and refreshing experienced emergency responders.

- Dynamic and versatile
- Easy-to-use and safe
- Bilingual training solution
- Same intuitive and simple 1-2-3-step operation as the Saver One Semi-Automatic AED
- Simulates real-life sudden cardiac arrest events with 10 preconfigured adult and pediatric scenarios activated by remote control
- Capable for training Advanced Pediatric CPR (15:2 CV ratio) in P-BLS advanced class
- LED control indicators for the real "low battery" status or for simulating a faulty unit event



Remote Control Buttons:

- Turn Off
- Pause/Play Scenario
- High/Low Volume
- Swap languages
- 10 Preconfigured Scenarios
- 6 Single Event Simulation
- 4 Customizable Scenarios (option)



Saver One T AED Trainer

Code: SVT-B0008

Comes with nylon carrying case, 2 pairs of training pads, 1 pads cable, 1 remote control, 1 rechargeable battery, 1 smart charger, 1 user's manual.

Dimensions: cm. 21,5 W x 26,5 D x 7,5 H

Weight: kg. 1,10 with battery

Battery: Ni-MH Rechargeable Battery (12V@2Ah) 100hrs

Charger: Smart Charger mod.H01440010-XX-D1 Ni-MH/Ni-Cd

Remote Control: powered by 2 AA 1.5V batteries (not provided)

1 Year Guarantee



info@amiitalia.com
www.amiitalia.com

Office & Production:

Via Cupa Reginella, 17/A
80010 Quarto
Tel: +39.081.8060574
Fax: +39.081.8764769

Office:

Viale Gran Sasso, 11
20131 Milano
Tel: +39.02.20509246
Fax: +39.02.29520839