

Baby Anne

This low-cost, lifelike infant CPR trainer is the perfect supplement to the Resusci Baby manikin

The Baby Anne manikin was developed to provide effective infant CPR training without compromising realism or quality. The convenient packaging of four Baby Anne manikins makes hands-on practice affordable for every student.

Product benefits:

- Educationally effective by offering all the essential features necessary for learning infant CPR
- Affordability allows increased hands-on practice
- Realistic chest compressions and chest rise allow students to learn proper techniques
- Logistically convenient lightweight and portable 4-pack weighs under 5,5 kg making Baby Anne easy to transport
- Practical because it is inexpensive and easy to maintain
- Durable construction allows unequalled long-term use

Product features:

- Natural obstruction of the airway allows students to learn the important technique of opening the airway
- Head tilt/chin lift and jaw thrust allow students to practice correctly all maneuvers necessary when resuscitating a real patient
- Realistic chest compliance means students can experience the proper technique required for chest compressions on infants
- Foreign-body airway obstruction feature allows the release of a foreign-body obstruction to be practiced through back blows and chest-thrust techniques
- Economical disposable airways for quick and easy clean up
- Removable and reusable faces for convenient and affordable maintenance



05 00 00 Baby Anne

05 00 02 Baby Anne Dark Skin

05 00 10 Baby Anne 4-Pack

05 00 12 Baby Anne 4-Pack Dark Skin



Release of foreign-body obstruction can be practiced.

Baby Anne Manikin Includes:

1 Manikin, Soft Pack, 6 Airways, 10 Foreign Body Practice Objects and Directions for Use.

Baby Anne 4-Pack Manikin Includes:

4 Manikins, 4-Pack Carry Case, 24 Airways, 40 Foreign Body Practice Objects and Directions for Use.

Consumables



05 02 00
Baby Anne Faces
(pkg. 6)



05 02 02
Baby Anne Faces,
dark skin (pkg. 6)



05 01 00
Airways (pkg. 24)



05 03 00
Foreign Objects
(pkg. 10)



15 12 01
Manikin Face Shields
(6 rolls of 36 each)